

Feature:

- 15 inch LED screen
- Pretty, thin, and light
- 2 probe connector
- Full digital imaging technology, crystal-clear Image
- 3D Image function
- PC based, abundant functions
- Detail report edit function
- Powerful image management
- Large volume storage image and cineloop
- Not only can save image but also can save video
- DICOM display and transfer
- Compatible with laser/inkjet printers
- Abundant calculation functions and can select different calculate formula
- Can save and select some groups inspect parameter so save time to adjust per time
- Built-in battery(Work 5 hours)
- Multilingual(English, French, German, Spanish and so on)
- Good aluminum trunk packing
- Multiple language: English, French, German, Spanish, Russian, Italian, Portuguese, etc.**

Specifications:

- Maximum Scanning Depth: max. 250mm
- Multi-frequency: 5 segment frequencies
- Cine Loop: max. 1024 frames
- Scanning mode: B, 2B, B/M, M, 4B, 9B, ZOOM;
- Real-time Zoom on B mode
- Image gray scale: 256 level gray scale
- Image processing: Pre-processing, after-processing, dynamic range, frame rate, line average, edge; enhancement, Black/White inversion; Gray scale adjustment, contrast, brightness, γ revision.
- Gain: 0-100dB; Time gain control(TGC): 8 segment adjustment, B, M adjustment separately.
- Measurement and calculation(human):
 - B mode: distance, circumference, area, volume, angle, ratio, stenosis, profile, histogram;
 - M mode: heart rate, time, distance, slope and stenosis;
 - Gynecology measurement: Uterus, cervix, endometrium, L/R ovary;
 - Obstetric: gestation age, fetal weight, AFI;
 - Cardiology: LV, LV function, LVPW, RVAWT;
 - Urology: transition zone volume, bladder volume, RUV, prostate, kidney;

Small parts: optic, thyroid, jaw and face.

Veterinary: Cattle: BPD, CRL, TD; Horse: GA; Dog: GS; Cat: HD, BD; Sheep: BPD, CRL, TD

-Supported output apparatus: U-Disk, **Laser/ink printer**, video printer, VGA

-Zoom: 10 ratio, 1.5, 2.0, 2.5, 3.0, 3.5, 4.0, 4.5, 5.0, 5.5, 6.0

-Image store: by hardisk or USB disk, **can save image and also video**(avi or swf format)

-**3D image function**(function): can get rebuilt 3D image by any probe

-Work time by built-in battery: about **5 hours**

-Monitor: **15 inch color LED screen**

-Probe Connector:2

-Dimension: 362mm×340mm×60mm

-Weight: Around 4.8Kg

Standard Configuration:

One main host, one convex probe(3.5MHz), one power adapter, one aluminum trunk.

Optional:

Linear, Trans-vaginal, Micro-convex,Rectal Probe; Bag; Trolley; Laser Printer.